



Create a project




LabVIEW version:

LabVIEW Full Development System
Copyright (c) 1986-2022 National Instruments. All rights reserved.
Version 2022 Q3 (64-bit) 22.3.4f8

Create a compilation



Requirement



Open the Package manager of labViewruntime et install the NI-VISA driver

Warning !

You have to set the COM port number for trigger communication. This is because this part of the code is outside the main loop, which means that the value at runtime is used. If you add a drop-down menu with the list of available COM ports, the value will simply be ignored because it is empty at the start. You must therefore set it (for example COM8).

